SAA-VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of SAA with solo analyser.

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Kit size

6 tests/pack
1 x 2 ml (decision level)

Test kit preparation: Allow single test min. 10 minutes to warm up to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.

Summary

Serum Amyloid A (SAA) is an early and sensitive blood biomarker for tissue injury and inflammation and has been indicated in many inflammatory diseases. The level of SAA circulating in the blood is known to increase dramatically in response to tissue damage or inflammation, classifying it as an acute phase protein. Circulating SAA concentrations may increase up to 1000-fold following inflammation, infection, tissue injury and cell necrosis and decline rapidly following recovery.

It has also been observed in several veterinary species including the horse that SAA is a very useful inflammatory marker that may be used for detection of clinical and possibly subclinical disease, monitoring of disease activity and response to therapy, prognostication, and detection of spread of infections in herds (Nielsen et al., 2004). The kinetic profile of the SAA response makes SAA an excellent indicator of inflammation. Hepatic SAA synthesis begins shortly after an inflammatory insult. As a result plasma concentrations start to increase within few hours, and they peak 36–48 h after injury, as demonstrated previously in horses after experimental induction of inflammatory or infectious disease (Hultén et al., 2002).

Method / Measurement Range

Latex agglutination test based on optical measurement of the change in turbidity caused by the agglutination of the latex particles sensitised with SAA antibodies.

Measurement ranges:

Horse: 10 - 500 µg/ml (10 - 500 mg/l)

Cat: 10 - 150 µg/ml (10 - 150 mg/l)

Samples with concentrations higher than the upper limit of the measurement range must be diluted 1 to 2 with physiological saline (0,9% NaCl), eg. 20 µl sample + 40 µl 0,9%-NaCl-solution, and the result then multiplied by 3.

Lower detection limit: 10 µg/ml

Stability and Storage

Store at 2 - 8°C. DO NOT FREEZE. Kit can be used until expiry date printed on the labels.

Sample Material

Use serum for routine analysis only. In several cases milk and synovial fluids can be used for SAA detection as well.

Warnings and Precautions

This testkit is for **veterinary** in vitro diagnostic use only! DO NOT INGEST. Avoid contact with skin and eyes. Contains sodium azide, which may react with lead or copper plumbin to form explosive compounds. Take the necessary precautions for the use of laboratory reagents.

Quality Control

Good laboratory practice suggests the use of control specimens to ensure proper assay performance. Controls results must fall within the established ranges, as determined by your laboratory. If results fall outside of the established ranges, assay results are invalid.

Reference Ranges

Cat: < 10 µg/ml

Horse: < 20 µg/ml

Camel: < 10 µg/ml

For further species please establish laboratory internal normal ranges. These ranges are given for orientation only. Each laboratory should establish its own reference values.

Precision for Plasma

cv (99 µg/ml) < 8 % n=20

Correlation

y = 0,98x + 0,12

y = "solo", x = (Eiken Chemicals SAA on Eurolyser CCA180) Reference Methode

Waste Management

Please refer to local legal requirements.

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English

ATTENTION!

Allow test min. 10 minutes to warm up to room temperature!


1. Preparation of test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette in test kit rack
- 1.3 Place ERS cap in test kit rack
- 1.4 Press measurement button, enter required information using the touch screen

2. For further details please see user manual of laboratory photometer

- 2.1 Start of analysis by pressing the start button.

3. Sample preparation serum or plasma

- 3.1a Pipette 5 µl serum from primary tube
OR ...
- 3.1b Pipette 5 µl plasma from centrifuged sample tube
- 3.2 Pipette 5 µl sample into ERS cuvette INTO LIQUID
- 3.3 Firmly apply ERS cap on ERS cuvette
- 3.4 Place ERS cartridge into laboratory photometer
- 3.5 Start automatic sample processing by pressing the  start button.